

```
In [1]: # lime_explain.py
# sudo -E /opt/tljh/user/bin/pip3 install lime

# conda activate jupyterlab-debugger38
```

```
In [1]: from lime.lime_image import *
import pandas as pd
import yaml #pyyaml
import os
import datetime
import dill
import cv2 #opencv-python
import numpy as np
from tensorflow.keras.models import load_model
from tensorflow.keras.preprocessing.image import ImageDataGenerator
import sys
#sys.path.append('/home/ubuntu/covid-cxr/src/')
sys.path.insert(0, os.path.abspath('../'))
from visualization.visualize import visualize_explanation
from predict import predict_instance, predict_and_explain
from data.preprocess import remove_text
```

```
In [2]: def setup_lime():
    ...
    Load relevant information and create a LIME Explainer
    :return: dict containing important information and objects for explanation experiments
    ...

    # Load relevant constants from project config file
    cfg = yaml.full_load(open("/home/ubuntu/covid-cxr/config.yml", 'r'))
    lime_dict = {}
    lime_dict['NUM_SAMPLES'] = cfg['LIME']['NUM_SAMPLES']
    lime_dict['NUM_FEATURES'] = cfg['LIME']['NUM_FEATURES']
    lime_dict['IMG_PATH'] = cfg['PATHS']['IMAGES']
    lime_dict['RAW_DATA_PATH'] = cfg['PATHS']['RAW_DATA']
    lime_dict['IMG_DIM'] = cfg['DATA']['IMG_DIM']
    lime_dict['PRED_THRESHOLD'] = cfg['PREDICTION']['THRESHOLD']
    lime_dict['CLASSES'] = cfg['DATA']['CLASSES']
    lime_dict['CLASS_MODE'] = cfg['TRAIN']['CLASS_MODE']
    lime_dict['COVID_ONLY'] = cfg['LIME']['COVID_ONLY']
    KERNEL_WIDTH = cfg['LIME']['KERNEL_WIDTH']
    FEATURE_SELECTION = cfg['LIME']['FEATURE_SELECTION']

    # Load train and test sets
    lime_dict['TRAIN_SET'] = pd.read_csv(cfg['PATHS']['TRAIN_SET'])
    lime_dict['TEST_SET'] = pd.read_csv(cfg['PATHS']['TEST_SET'])
```

```

# Create ImageDataGenerator for test set
test_img_gen = ImageDataGenerator(preprocessing_function=remove_text,
                                   samplewise_std_normalization=True,
                                   samplewise_center=True)

test_generator = test_img_gen.flow_from_dataframe(dataframe=lime_dict['TEST_SET'],
directory=cfg['PATHS']['RAW_DATA'],
          x_col="filename", y_col='label_str', target_size=tuple(cfg['DATA']['IMG_DIM']),
batch_size=1,
          class_mode='categorical', validate_filenames=False, shuffle=False)
lime_dict['TEST_GENERATOR'] = test_generator

# Define the LIME explainer
lime_dict['EXPLAINER'] = LimeImageExplainer(kernel_width=KERNEL_WIDTH,
feature_selection=FEATURE_SELECTION,
                                   verbose=True)

dill.dump(lime_dict['EXPLAINER'], open(cfg['PATHS']['LIME_EXPLAINER'], 'wb')) #
Serialize the explainer

# Load trained model's weights
lime_dict['MODEL'] = load_model(cfg['PATHS']['MODEL_TO_LOAD'], compile=False)

return lime_dict

```

In [5]:

```

def explain_xray(lime_dict, idx, save_exp=True):
    ...

    Make a prediction and provide a LIME explanation
    :param lime_dict: dict containing important information and objects for explanation
experiments
    :param idx: index of image in test set to explain
    :param save_exp: Boolean indicating whether to save the explanation visualization
    ...

    # Get i'th preprocessed image in test set
lime_dict['TEST_GENERATOR'].reset()
for i in range(idx + 1):
    x, y = lime_dict['TEST_GENERATOR'].next()
x = np.squeeze(x, axis=0)
x

# Get the corresponding original image (no preprocessing)
orig_img = cv2.imread(lime_dict['RAW_DATA_PATH'] + lime_dict['TEST_SET']['filename'][idx])
new_dim = tuple(lime_dict['IMG_DIM'])
orig_img = cv2.resize(orig_img, new_dim, interpolation=cv2.INTER_NEAREST) # Resize
image

# Make a prediction for this image and retrieve a LIME explanation for the prediction
start_time = datetime.datetime.now()
#explanation = explainer.explain_instance(images[0].astype('double')) # added double for

```

```

tensorflow2_latest_p37
    explanation, probs = predict_and_explain(x.astype('double'), lime_dict['MODEL'],
lime_dict['EXPLAINER'],

                                                lime_dict['NUM_FEATURES'],
lime_dict['NUM_SAMPLES'])
    print("Explanation time = " + str((datetime.datetime.now() - start_time).total_seconds())
+ " seconds")

    # Get image filename and label
    img_filename = lime_dict['TEST_SET']['filename'][idx]
    label = lime_dict['TEST_SET']['label'][idx]

    # Rearrange prediction probability vector to reflect original ordering of classes in
project config
    probs = [probs[0][lime_dict['CLASSES'].index(c)] for c in
lime_dict['TEST_GENERATOR'].class_indices]

    # Visualize the LIME explanation and optionally save it to disk
    if save_exp:
        file_path = lime_dict['IMG_PATH']
    else:
        file_path = None
    if lime_dict['COVID_ONLY'] == True:
        label_to_see = lime_dict['TEST_GENERATOR'].class_indices['COVID-19']
    else:
        label_to_see = 'top'
    _ = visualize_explanation(orig_img, explanation, img_filename, label, probs,
lime_dict['CLASSES'], label_to_see=label_to_see,
                            dir_path=file_path)

    return

```

```

In [8]: if __name__ == '__main__':
        lime_dict = setup_lime()
        i = 3 # Select i'th image in test set
        explain_xray(lime_dict, i, save_exp=True) # Generate explanation for
image

```

Found 182 non-validated image filenames belonging to 2 classes.

```



Intercept -0.03701115675992887
Prediction_local [0.19896746]
Right: 0.22185671
Intercept 1.0370111510015623
Prediction_local [0.80103255]
Right: 0.7781433
Explanation time = 8.989423 seconds

```

Ground Truth Class: 0 (non-COVID-19)

Predicted Class: 0 (non-COVID-19)

Prediction probabilities: ['0.78', '0.22']



In [ ]: